

Industrial Pneumatic Tube Transport System

EFFICIENT SAMPLE TRANSPORT

Our pneumatic airtube systems ensure the reliable and swift transportation of samples from the steelworks to the laboratory. Leveraging our extensive decades-long expertise in this field, we are adept at installing these systems even in challenging spatial conditions and across vast distances.

Sample data transfer

Utilizing our advanced sending and receiving stations, we seamlessly integrate sample information from the customer's manufacturing execution systems (MES) into our Sample Control Software (SamCoS). Sample registration can be operated easily by using buttons our touch-display or RFIDtechnology (by scanning a chip within the carrier).

Manual or fully automated

Upon the pneumatic tube's arrival at the laboratory, the sample undergoes a (1) manual or automatic unpacking process and is carefully aligned for subsequent sample preparation and analysis. Powdered samples, such as slag, are discharged into cups for further processing. This step occurs either within the (2) receiving station itself or at (3) specialized unpacking and transfer stations. In every scenario, our processes are meticulously optimized to achieve short throughput times and ensure maximum plant availability.



1Td 1 VSdee1 1 Ve1I S@cQdeI 1eT1Vc1







Manual Sending Station SAMPLI-SRM-M

RELIABLE AND FAST SENDING

The SAMPLI-SRM-M serves as our plant station designed for the seamless sending and receiving of airtube carriers, facilitating the transit from the steelworks to the laboratory. The steel, iron, or slag sample is carefully loaded into the designated carrier, sealed with the aid of the integrated tool, and positioned at the input point. Once the sample identification has been chosen through a simple pushbutton on the control panel, the carrier embarks on its journey to the laboratory, where it can be unpacked for further processing. Optionally automated sample ID recognition is possible by RFID technology. The station is sesigned for high sample throughput and 24/7 shift operation. Our design meets safe and simple operation combined with special safety measures. It is characterized by particular robustness and reliability, especially in the harsh environment of production plants.

Technical specifications

PRODUCT DATASHEET

USE	Sending and Receiving station for carriers within production facility		
Drive		3x 400V, 50 HZ	
Protec	Protection Class IP 54		
Dimensions (WxDxH): Approx. 960x400		Approx. 960x400x1.915mm	
Weight		Approx. 140 kg	





1Td 1 VSdee1 1 Ve11 S@Qdel 1eT1Vc1 Carrier-Opener

Operator Terminal with buttons



Auto Receiving Station SAMPLI-RSM-UNPK

FLEXIBLE AND EFFICIENT HANDLING WITHIN STATION

The SAMPLI-RSM-UNPK dispatch and receiving station is specifically engineered for the automated handling of pneumatic tube carriers within the laboratory context. This station facilitates the receiving and sending of diverse sample materials, such as hot or cold metal samples, granular materials, powder samples, and other types of samples. Upon receiving a pneumatic tube carriers, the this system automatically performs a sequence of actions, including opening, emptying, closing, and returning the carriers to the dispatching station. This entire process is seamlessly managed by an integrated handling system. All system modules are seamlessly integrated into the machine column and are easily accessible through doors in the housing.

Technical specifications

USE	Opening,	emptying	within	laboratory	and	closing,	and
	returning	the carriers	to the f	ield station.			
Drive	Drive 3x 400V, 50 HZ						
Protect	Protection Class IP 54						
Dimensions (WxDxH): Approx. 770x500x2.23		ōmm					
Weight			Approx.	290 k	g		





Integrated handling system



1Td 1 VSdee1 1 Ve1I S@cQdel 1eT1Vc1



Auto Receiving Station SAMPLI-RSM-R

ROBUST & COMPACT DESIGN FOR DIRECT WAYS TO LAB

The SAMPLI-RSM-R laboratory station is significantly more space-saving than the SAMPLI-RSM-UNPK in order to enable as many direct routes (transport lines without diverters) to the laboratory as possible. SAMPLI-RSM-R is specifically designed for the automated dispatch and reception of pneumatic carriers in robot automation. All integrated machines are arranged in a circular fashion and are operated by a central robot. Upon the carrier's arrival in the laboratory, the robot extracts it from this station and transfers it to the unpacking station SAMPI-OPST-R (see page 6). Here, the sample is removed for further processing in the automation system. Subsequently, the robot returns the empty carrier to SAMPLI-RSM-R to send it back to the plant station. If necessary, operators can use this station for manual operation from back site.

Technical specifications

USE	Receipt of pneumatic carriers in the robot laboratory and		
	sending back to the field station		
Drive 3x 400V, 50 HZ		3x 400V, 50 HZ	
Protec	ction Class	IP 54	
Dimensions (WxDxH): Approx. 430x200x1.		Approx. 430x200x1.900mm	



Front site

Back site (for manual operation)



1Td 1 /Sdee1 1 Ve1I SQeQdel 1eT1Vc1



Tube Diverter and Tube Transport Contact



RELIABLE AND FAST TRANSPORT

DIVERTER redirect carriers from two or four tubes along one selected route to the laboratory (e.g. From 2 or 4 SENDING STATIONS to 1 LABORATORY RECEIVING STATION). Conversely, they can be used to split a pneumatic tube line into two or four tubes. TUBE TRANSPORT CONTACT is used to signal Tube Transport Contact and control carrier that have passed through.

2-Way / 4-Way Diverter

PRODUCT DATASHEET

TRANSPORT BLOWER STATION



SUITABLE FOR LONG DISTANCES

The TRANSPORT BLOWER STATION generates driving air for the airtube system. Depending on the direction of rotation of the motor (for blowers with max. 3.6 kW) or the position of the reverse valve (for blowers up from 4.0 KW engine power), driving air is fed into the pneumatic tube system or withdrawn from it, thus determining the direction of travel of a sample carrier. This system shortens sampling times and increased efficiency due to return valves that allow rapid change of air direction over several kilometers (no need to change the direction of rotation of the electric motor). Moreover the robust design ensures reliable transport even under particularly harsh environments.

Technical specifications

Engine power	3,6kW, 4,8kW, 7,5kW or 11kW	
Weight	56kg, 90kg, 108kg or 150kg	
Dimensions (LxWxH): 484 x 595 x 583mm		
	750 x 600 x 868mm	
	750 x 600 x 993 mm	
	750 x 600 x 1.100m	



1Td 1 VSdee1 Ve1I SceQdel 1eT1Vc1

Carrier-Opener and Conveyor Belt

SAMPLI-OPST-R (STEEL) SAMPLI-OPSL-R (SLAG)



Technical specifications

Dimensions (WxDxH): Approx. 300x300x1.820mm		
Drive	Pneumatic	
USE	Opening of Sample Carrier Steel	



Technical specifications

USE	Opening of Sample Carrier Steel
Drive	Pneumatic
Dimensions (WxDxH): Approx. 300x300x1.615mm	

CONVEYOR BELT



HIGHLY CONFIGURABLE

This conveyor belt is specially configurable for transporting sample cups, steel rings or samples itself within the process chain of quality laborator This conists the transport to sample preparation machines as well as to XRF, XRD or OES analyis units. This system is fully automated through communication to higher-level control.

Technical specifications

Operating voltage	24V
Weight	Project dependent

Dimensions (WxDxH): Project dependent



1Td 1 VSdee1 1 Ve1I S@cQdeI 1eT1Vc1